

Detection of KMT2A-PTDs in healthy donor and myeloid malignant samples using next generation sequencing

Jiajie Huang^a, Haigang Gu^b, Janet Orton^b, Marina Sedova^a, Amir Marcovitz^a, Jennifer Burke^a, Sarah Brozio^a, Paul Williams^c, Scott Myrand^c, Nate Olowo^a, Adam Broomer^d, Brendan Deal^a, Collyn Seeger^e, Seth Sadis^c, Sophie Rozenhak^d, Fiona Hyland^a, Guang Liu^b
 a: Thermo Fisher ScientificTM, South San Francisco, CA; b: Sonora Quest Laboratories, Phoenix, AZ; c: Thermo Fisher ScientificTM, Ann Arbor, MI; d: Thermo Fisher ScientificTM, Carlsbad, CA; e: Thermo Fisher ScientificTM, Guilford, CT

INTRODUCTION

KMT2A (*MLL*) fusions and *KMT2A*-PTD (partial tandem duplication) are vital biomarkers in myeloid malignancies traditionally detected by RT-qPCR (quantitative real-time PCR). This study utilizes next generation sequencing (NGS) with OncoPrintTM Myeloid Assay GX v2 to report the detection of *KMT2A*-PTDs in both healthy donors and myeloid malignancy samples. *KMT2A* fusions in myeloid malignant samples are also reported.

Fig 1. Estimated new cases (%) of Leukemia, Lymphoma and Myeloma in USA, 2021




Fig 2. OncoPrintTM Myeloid assay enables rapid and efficient multi-biomarker testing




Fig 3. Genexus Instruments



MATERIALS AND METHODS

We sequenced 8483 research samples with known myeloid malignancies in both Sonora Quest LaboratoriesTM and at Thermo Fisher ScientificTM South San Francisco site. We also acquired 20 healthy donor whole blood samples (total 127 replicates) from StanfordTM Blood Center and Discovery Life SciencesTM and sequenced them at 3 different sites of Thermo Fisher ScientificTM (South San Francisco, CA; Guilford, CT; Carlsbad, CA). Samples were processed on the Ion TorrentTM GenexusTM Software 6.6 and analyzed using the OncoPrintTM Myeloid Assay GX v2 for fusion profiling targeting 6 different *KMT2A*-PTD variants and 199 *KMT2A* fusion isoforms.

Table 1. OncoPrintTM Myeloid Assay GX v2 Panel

DNA Panel		RNA Panel				
Hotspot genes (28)	Full genes (17)	Fusion Driver Genes (30)	Expression genes (5)	Expression control genes (5)		
ANKRD26 ABL1 BRAF CBL CSF3R DDX41 DNMT3A FLT3 (ITD + TKD) GATA2 HRAS IDH1 IDH2 JAK2 KIT	KRAS MPL MYD88 NPM1 NRAS PPM1D PTPN11 SMC1A SMC3 SETBP1 SF3B1 SRSF2 U2AF1 WT1	ASXL1 BCOR CALR CEBPA ETV6 EZH2 IKZF1 NF1 PHF6 PRPF8 RB1 RUNX1 SH2B3 STAG2 TET2 TP53 ZRSR2	ABL1 BCL2 BRAF CCND1 (MKL1) CREBBP EGFR FGFR1 FGFR2 HMGA2 JAK2 KAT6A (MOZ) KAT6B KMT2A (MLL) RARA RUNX1 TCF3 TFE3 ZNF384	MECOM MET MRTFA (MYBL1) MYBL1 NTRK2 NTRK3 NUP214 NUP98 PAX5 PDGFRA PDGFRB	BAALC MECOM MYC SMC1A WT1	EIF2B1 FBXW2 PSMB2 PUM1 TRIM27

Fig 4. *KMT2A*-PTDs: a relevant biomarker in myeloid malignancies



Table 2. OncoPrintTM Myeloid assay is designed to detect six *KMT2A*-PTD variants

Variant	Length (bp)
<i>KMT2A-KMT2A</i> .nt51.K11K2	178
<i>KMT2A-KMT2A</i> .K9K2	172
<i>KMT2A-KMT2A</i> .K7K2	169
<i>KMT2A-KMT2A</i> .K11K2	116
<i>KMT2A-KMT2A</i> .K8K2	105
<i>KMT2A-KMT2A</i> .K10K2	105

Table 3. Report only the PTDs associated with myeloid malignancy with high confidence

Example: *KMT2A*-PTDs @ Genexus software 6.6 OMAv2 w4.3.2 (a past release)

Sample name	Isoforms detected	Read count
Healthy donor sample #A	<i>KMT2A-KMT2A</i> .K8K2	9
	<i>KMT2A-KMT2A</i> .K10K2	147
	<i>KMT2A-KMT2A</i> .K7K2	201
	<i>KMT2A-KMT2A</i> .K9K2	343
Myeloid malignant sample #B	<i>KMT2A-KMT2A</i> .K10K2	551
	<i>KMT2A-KMT2A</i> .K7K2	722
	<i>KMT2A-KMT2A</i> .K8K2	311
	<i>KMT2A-KMT2A</i> .K9K2	2466

Is there a differentiation in PTD read count in myeloid malignant vs healthy donor samples?

RESULTS


The mean read length of this data set is 90 – 120 bp and the mapped fusion reads is 20,000 – 30,000. *KMT2A*-PTDs were detected in both healthy donors and myeloid samples. Healthy donor PTD read counts were consistently <2000 and averaged 1/3 of myeloid samples. About 33% of myeloid samples had higher PTD read counts than any healthy donor sample. BLAT (BLASTTM-like Alignment Tool) analysis confirmed specific exon matching on the *KMT2A* gene in both cohorts. Among the 8503 myeloid samples, 162 contained a total of 5 unique *KMT2A* PTDs, and 105 contained a total of 30 unique *KMT2A* fusion isoforms with *KMT2A-MLLT1* and *KMT2A-MLLT3* being the most prevalent *KMT2A* fusion gene pairs.

Table 4. *KMT2A* fusions & PTDs existed in ~3% samples

	Gene Pair	Found in # of unique samples	Found in % of all samples (N=8503)
<i>KMT2A</i> fusions	<i>KMT2A-AFF4</i>	2	0.02%
	<i>KMT2A-CASC5</i>	2	0.02%
	<i>KMT2A-CBL</i>	4	0.05%
	<i>KMT2A-ELL</i>	11	0.13%
	<i>KMT2A-EPS15</i>	4	0.05%
	<i>KMT2A-MLLT1</i>	30	0.35%
	<i>KMT2A-MLLT10</i>	9	0.11%
<i>KMT2A</i> -PTDs	<i>KMT2A-MLLT3</i>	29	0.34%
	<i>KMT2A-MLLT4</i>	14	0.16%
	Total	105	1.23%
	Total	162	1.91%
Samples w/ ≥ 1 <i>KMT2A</i> fusion or PTD	Total	265	3.12%

~2% samples had at least one *KMT2A*-PTD. ~1% samples had at least one *KMT2A* fusion.

Fig 5. Read count differentiates myeloid cancer PTDs



REFERENCES

Schnitter S, Wörmann B, Hiddemann W, Griesinger F. Partial tandem duplications of the *MLL* gene are detectable in peripheral blood and bone marrow of nearly all healthy donors. *Blood*. 1998 Sep 1;92(5):1728-34. PMID: 9716602.
 Dai B, Yu H, Ma T, Lei Y, Wang J, Zhang Y, Lu J, Yan H, Jiang L, Chen B. The Application of Targeted RNA Sequencing for *KMT2A*-Partial Tandem Duplication Identification and Integrated Analysis of Molecular Characterization in Acute Myeloid Leukemia. *J Mol Diagn*. 2021 Nov;23(11):1478-1490. doi: 10.1016/j.jmoldx.2021.07.019. Epub 2021 Aug 10. PMID: 34384895.

Fig 5. IGV view of *KMT2A-KMT2A.K9K2* in healthy donor samples vs myeloid samples



On IGV view (Fig 5(a)), the alignment of *KMT2A-KMT2A.K9K2* looks clean with some insertions and mismatches in both healthy donor and myeloid malignant samples. This is confirmed by BLAT (Fig 5(b)). One insertion of 18 bps is observed at the break point of *KMT2A-KMT2A.K9K2*. BLAT shows that it is from intron 1 (Fig 5(c)).

CONCLUSIONS

In this study, we describe the detection of *KMT2A* fusions in myeloid malignant samples. Our study also describes the detection of *KMT2A*-PTDs in both healthy donor and myeloid samples, with myeloid cases showing significantly higher PTD read counts. Additional studies to understand the relevant expression level of PTD are in progress. This intriguing finding opens opportunities for prospective studies to monitor individuals with elevated PTD levels for myeloid malignancy development and retrospective studies to explore whether healthy donors identified with this alteration years ago after blood donation were subsequently recorded in the national health system with myeloid malignancies.

ACKNOWLEDGEMENTS

Thank you to all the individuals involved in the development of the OncoPrintTM Myeloid Assays.

DISCLAIMER

For Research Use Only. Not for use in diagnostic procedures. © 2024 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Sonora Quest Laboratories is a trademark of Sonora Quest Laboratories, LLC. Discovery Life Sciences is a trademark of Discovery Life Sciences, LLC. BLAST is a trademark of The National Library of Medicine.

