

# Increasing Positive Clonality Detection Rate in Multiple Myeloma Research Samples Using NGS Characterization of Multiple B cell Receptors in a Single Reaction


Geoffrey M. Lowman<sup>1</sup>, Landon Pastushok<sup>3</sup>, Karen Mochoruk<sup>3</sup>, Wayne Hill<sup>3</sup>, Michelle Toro<sup>1</sup>, Shrutii Sarda<sup>2</sup>, Chenchen Yang<sup>2</sup>, Julie Stakiw<sup>4</sup>, Mark Bosch<sup>4</sup>, Hadi Goubran<sup>4</sup>, Ron Geyer<sup>3</sup>, John DeCoteau<sup>3</sup>


- [1] Clinical Sequencing Division, Thermo Fisher Scientific, Carlsbad, CA 92008, USA [2] Clinical Sequencing Division, Thermo Fisher Scientific, South San Francisco, CA 94080, USA
- [3] Advanced Diagnostics Research Laboratory, University of Saskatchewan, Health Sciences 4D01.12, Saskatoon SK, Canada [4] Saskatoon Cancer Agency, 20 Campus Dr, Saskatoon, SK S7N 4H4, Canada

#### **INTRODUCTION**

B cell repertoire analysis by next-generation sequencing (NGS) is at the forefront of leukemia and lymphoma research. Some advantages provided by NGS-based techniques include a lower limit-of-detection and simpler paths to standardization compared to other methods. Importantly, in research of post-germinal B cell disorders, such as multiple myeloma, NGS methods allow for the study of clonal lineage based on somatic hypermutation (SHM) patterns. Current targeted NGS assays require multiple libraries to survey each B cell receptor chain (IGH, IGK, IGL), and this fact is highlighted when initial clonality detection fails due to mutations under primer binding sites. This issue can be especially true with multiple myeloma which has a high rate of SHM. To address these issues, we have developed an assay for B cell analysis, based on Ion AmpliSeq™ technology, which enables efficient detection of IGH, IGK, and IGL chain rearrangements in a single reaction.

## **PANEL DESIGN**





FR2-J

CDR3

~250bp

Diversity (D)

CDR1

Secondary

Assays

CDR2

#### **MATERIALS AND METHODS**

**Commercially procured** 

The BCR Pan-Clonality panel targets the framework 3 (FR3) portion of the variable gene and the joining gene region of heavy- and lightchain loci (IGH, IGK, IGL) for all alleles found within the IMGT database, enabling readout of the complementary-determining region 3 (CDR3) sequence of each immunoglobulin chain. To maximize sensitivity, we included primers to amplify IGK loci rearrangements involving Kappa deletion element and the constant region intron. To evaluate assay performance, we conducted reproducibility studies and clonality assessment using gDNA from a total of 45 multiple myeloma research samples. All multiple myeloma cases examined in this work were confirmed clonal previously by light chain restriction via flow cytometry or IHC/ISH in tissue sections. Sequencing and clonality analysis was performed using the Ion GeneStudio S5 System and Ion Reporter 5.16 analysis software.

#### **Multiple Myeloma cohort:** Goal: Assign clonotype to respective MM 45 Clinical research BMA specimens research samples through detection of clonal BCR rearrangements Clinically characterized Wide range of disease burden by Following clonality assignment guidelines from: flow cytometry and/or IHC Arcila, Maria E. et.al. Establishment of Immunoglobulin Heavy (IGH) Chain Clonality 16 lambda light chain Testing by Next-Generation Sequencing for positive samples Routine Characterization of B-Cell and Plasma Cell Neoplasms. Journal of Molecular Diagnostics, Vol. 21, No. 2, March 2019. Included known polyclonal control Achieve clonal rearrangement detection in >85% MM samples

## **RESULTS**

## Oncomine<sup>™</sup> BCR Pan-Clonality Assay

Samples	Total Tested	Total positive (IGH) [%]	Total positive (IGL) [%]	Total positive (IGH+IGL) [%]	
		[ 70]	[ 70]	[ 70]	
Clinical Research Samples (MM)	45	34 [76%]	41 [91%]	42 [93%]	

### Oncomine<sup>™</sup> IGH FR3(d)-J and FR2-J Assays

Samples	Total Tested	Total positive FR3(d)-J [%]	Total positive FR2-J [%]
Clinical Research	14**	7	7
Samples (MM)		[50%]	[50%]

<sup>\*\*</sup> Testing included 11 samples not detected with the Pan-Clonality (IGH) and 3 borderline cases

#### Comparison of results from our MM (Post-GC) cohort (N=45) with a cohort of B-ALL (Pre-GC) samples (N=11)

## Both cohorts showed potential biased IGHV usage with V3-11 in 5/45 MM samples & Stereotypy IGHV4-34 was detected in 5/45 MM cases (and one B-ALL case) The high prevalence of IGHV4-34 in several lymphoma subtypes (e.g. BL, CLL, MZL) has led to a suggestion that autoreactive antigens may drive the expansion of some B-cell lymphomas.

SHM\*

All but two MM samples have SHM while all B-ALL samples have 0% SHM. When sequenced with FR1-J primers - low SHM in two MM samples (0.44% and 1.3%)

Most MM samples have high levels of SHM (6/35 > 10% SHM) Ion Reporter lineage analysis tool we identified 8/45 MM cases with 5 or more clones - one case with 23 clones!

In two samples only a clonal lambda rearrangement was detected.

clonal samples

14/16 (88%)

lambda-positive MM samples

were deemed clonal by the

IgK/IgL assay component

A clonal lambda light chain

rearrangement was identified

in 10/14 (71%) of the lgK/lgL

# Frequency of commercially available control detected by BCR Pan-Clonality Assay

	Polyclonal sample	Clonal sample 1	Clonal sample 2	Clonal sample 3
BCR Pan				
Clonality (lgK)	0.79 x 10 <sup>-4</sup>	0.78 x 10 <sup>-4</sup>	0.83 x 10 <sup>-4</sup>	0.87 x 10 <sup>-4</sup>
(19.1)				
BCR Pan- Clonality (IgH)	1.14 x 10 <sup>-4</sup>	3.47 x 10 <sup>-4</sup>	4.88 x 10 <sup>-4</sup>	2.35 x 10 <sup>-4</sup>
	Clonality (IgK) BCR Pan- Clonality	BCR Pan Clonality (IgK)  BCR Pan- Clonality 1.14 x 10 <sup>-4</sup>	BCR Pan Clonality (IgK)  BCR Pan-Clonality  1.14 x 10 <sup>-4</sup> 3.47 x 10 <sup>-4</sup>	BCR Pan Clonality (IgK)  BCR Pan-Clonality  1.14 x 10 <sup>-4</sup> 1.15 x 10 <sup>-4</sup> 1.16 x 10 <sup>-4</sup> 1.17 x 10 <sup>-4</sup> 1.18 x 10 <sup>-4</sup>

## **CONCLUSION**

Clonality assessment of multiple myeloma clinical research samples show a 93% overall positive detection rate by an assay which combines the IGH, IgK, and IgL chains in a single reaction using published guidelines for clonality assignment. 34 of 45 samples show positive detection of an IGH rearrangement, while 41 of 45 showed positive detection of at least one light chain receptor. In total, 42 of 45 samples were deemed clonal by the single tube assay based on detection for one or more receptor. Clonality results for this sample set are well correlated with orthogonal data from flow, IHC/ISH, or alternate NGS technologies. These results demonstrate the utility of a novel Ion AmpliSeq-assay for combined clonality analysis of B cell receptor heavy and light chains. We expect this assay to simplify workflow and open new paths for research in B cell disorders.

### **CORRESPONDENCE**

**Geoffrey Lowman** geoffrey.lowman@thermofisher.com



<sup>\*</sup>Oncomine BCR Pan-Clonality reports somatic hypermutation rate for the portion of the IGHV gene covered by the assay.