




# Oncomine Precision Assay on Ion Torrent Genexus System

Andy Felton

The Genexus System and its assays are for Research Use Only. Thermo Fisher Scientific is not promoting or encouraging any current diagnostic or clinical use of the Genexus System or any other Research Use Only products discussed in this presentation





# Cancer is a Disease of the Genome Caused by its Alterations



An effective and efficient tool is required to interrogate the alterations that cause cancer

\* <http://humanbiologylab.pbworks.com/w/page/48192744/Genomes%20%20Finding%20a%20Cure%20for%20Cancer>


# Biomarker Development is Accelerating



73% of medicines in oncology pipelines have associated biomarkers\*

# NGS is a Foundation of Precision Oncology Clinical Research

NGS can detect many different types of biomarkers simultaneously from a single sample



# Today's Challenges and Barriers for NGS Implementation in a Broader Lab Spectrum

## Too slow




Requires days and often weeks to get the results

## Too complex



High level of user expertise required to run NGS  
Modular workflows requiring multiple instruments and touchpoints

## Too costly



Cost of hiring and training staff  
Cost penalty for running small sample batches

## Too limited



Tissue requirements / QNS (quantity not sufficient) related failures



# A New Day for Genomic Profiling




**ThermoFisher**  
SCIENTIFIC

The world leader in serving science

# Oncomine Precision Assay on Ion Torrent Genexus System

A new generation solution for genomic profiling

## Fast



Single day  
sample-to-report.

## Hands free




2 touch points and only 10  
min of hands-on time.

## Cost saving



No more need to batch  
samples across multiple  
sequencing runs.

## Tissue saving



Minimum sample  
requirement  
Maximum results obtained.

\*Specimen-to-report workflow will be available after the Ion Torrent™ Genexus™ Purification System and integrated reporting capabilities are added in 2020.  
The content provided herein may relate to products that have not been fully validated by Thermo Fisher Scientific and is subject to change without notice.

# Introducing Ion Torrent™ Genexus™ System

The world's first turnkey, automated NGS system

Specimen to report in a single day with only two user touchpoints

- ✓ **Single-day turnaround time** allowing you to provide IHC and NGS results at the same time
- ✓ **Automated** sample purification, library prep, sequencing, and analysis reporting helps increase reproducibility and efficiency, while reducing personnel costs
- ✓ **Flexibility of economically running few or one sample** reduces the need for batching and helps you to deliver results faster than ever before
- ✓ System and consumables manufactured at a **facilities registered with FDA and ISO 13485 certified**
- ✓ **Research applications include** oncology & inherited, with reproductive health and infectious to come in the future



\*Specimen-to-report workflow will be available after the Ion Torrent™ Genexus™ Purification System and integrated reporting capabilities are added in 2020. The content provided herein may relate to products that have not been fully validated by Thermo Fisher Scientific and is subject to change without notice.




# Genexus System—Tomorrow's Specimen-to-Report NGS Workflow

## Genexus Software

Nucleic acid purification  
and quantitation\*

**Ion Torrent™ Genexus™**  
Purification System (Available 2020)



2 hour turnaround time  
12 FFPE (DNA and RNA)  
6 Plasma

Library preparation to  
variant interpretation

Report\*

**Ion Torrent™ Genexus™**  
Integrated Sequencer (Available November 2019)

**Ion Torrent™  
GX5™ Chip:**  
12–15M  
reads/lane



14 hours for a single-lane run  
(approx. 24 to 30 hours for full chip)  
Up to 32 Samples per run

- FFPE tissue
- Frozen tissue
- Bone marrow
- Whole blood
- PBL
- Urine
- Saliva

\*Specimen-to-report workflow will be available after the Ion Torrent™ Genexus™ Purification System and integrated reporting capabilities are added in 2020. The content provided herein may relate to products that have not been fully validated by Thermo Fisher Scientific and is subject to change without notice.

# Genexus Software—End-To-End Integration from Specimen to Report



## Integrated

Fully integrated solution enabling specimen-to-report workflow; no Ion Reporter server required



## Easy to use

Simplified, new user experience helps minimize the learning curve and human error



## Robust

Benchmarks on variant calling accuracy



## Flexible

Option to choose between integrated analysis on instrument or analysis on IR server or cloud

**example Labs**

Example Labs  
123 Street  
City, ST 12345 USA  
Tel (123) 123-1234  
email@example.com  
www.example.com

Sample ID: 00-123456789 Date: 05 Nov 2019 1 of 6

**Sample Information**

Year of Birth:	1986	Primary Tumor Site:	Stomach
Gender:	Male	Sample Type:	FFPE
Smoking Status:	Smoker	Sample ID:	1234
Case Number:	00-1234	Date Collected:	02-01-18

**Sample Cancer Type:** Gastric Cancer

**Variant Findings** ■ Indicated ■ Contraindicated

Genomic Alteration	Relevant Therapies (in this cancer type)	Relevant Therapies (in other cancer type)	Clinical Trials
FGFR1 fusion Tier: IC	None	None	5

Sources Included in relevant therapies: FDA<sup>1</sup>, NCCN, EMA<sup>2</sup>, EMA<sup>3</sup>

Disclaimer: The data presented here is a result of the curation of published data sources, but may not be exhaustive. The data version is 2019.09.000.

\*Specimen-to-report workflow will be available after the Ion Torrent™ Genexus™ Purification System and integrated reporting capabilities are added in 2020.  
The content provided herein may relate to products that have not been fully validated by Thermo Fisher Scientific and is subject to change without notice.

# Genexus™ Integrated Sequencer Maximizes Flexibility


- ✓ **Up to four different assays** prepared and sequenced simultaneously in a single run
- ✓ **Multiplexing capability** of up to 32 library preparations in a single run  
32 single-pool | 16 two-pool | 8 four-pool
- ✓ **Two week** on-instrument chip and reagent stability
- ✓ **Assay and throughput flexibility** facilitated by strips and multi-lane sequencing chip
- ✓ **Minimized consumable footprint** reduces required storage space





The content provided herein may relate to products that have not been fully validated by Thermo Fisher Scientific and is subject to change without notice.



## Library strips 1 and 2



## Templating strips 3 and 4



Strips 1 through 4, with descriptors for library chemistry or chip type, are installed in corresponding reagent bays on the deck of the Integrated Sequencer



**Ion GX5 Chip:**  
12–15M reads/lane



**Genexus™ Cartridge**

## Comparison of 4 Workflows

- DL8/Chef/S5
- Manual Library Prep/Chef/S5
- Manual Library Prep/ MiSeq
- Genexus

## Users

- 2 novice users and 3 expert AmpliSeq users

## Training


- All users received overview training on ILMN workflow
- All users are novice users of Genexus

## Technical Decisions

- 8 libraries were created for each workflow by each operator
- Decision not to use Biomek – customers will most likely not use Biomek for only 8 libraries





# Comparisons – Total Hands-On Time (Does not include Sample QC)



Users	Manual/Chef/S5	Manual/MiSeq	Genexus
Operator 1	01:50	02:05	N/A
Operator 2	03:01	03:49	00:12
Operator 1	01:12	02:34	00:12
Operator 2	02:24	03:26	00:14
Operator 3	01:22	02:34	00:14
Average	01:56	02:54	00:13

# Workflow Micro-costing Study Results




- Samples: 8 samples per run (OFA)
- FTE rate: \$33/hr (Glassdoor and Indeed rates for Clinical Lab Tech)
- Genexus FTE cost per run: \$7.15
- MiSeq FTE cost per run: \$93.61 (short of data analysis FTE costs)

Difference between Genexus and MiSeq FTE cost

- 52 weeks
- 4 year based on amortization of capital assets

# Genexus Is Cost Effective at Low Batch Sizes

Compared to GeneStudio, Genexus has steadier price per sample across a wide range of batch sizes which allows laboratories to test samples without having to always wait for a large batch size



#### Assumptions:

- For Genexus, per sample pricing calculated using list prices for Oncomine Precision Assay with associated Genexus consumables (running FFPE samples)
- For GeneStudio, per sample pricing calculated using list prices for Oncomine Focus Assay with associated GeneStudio S5 consumables (Ion Chef automated library preparation)



# Genexus Integrated Sequencer and Assay Turnaround Times

Metric	Assay	1 Lane	2 Lanes	3 Lanes	4 Lanes
Output in Reads (M)		12-15	24-30	36-45	48-60
Turnaround Time ( <i>Batch Size</i> )  Nucleic Acid to Oncomine Report	Custom AmpliSeq 1 Pool Assay	14 (4)	NA	NA	24 (16)
	Oncomine Comprehensive Assay v3 (DNA & RNA workflow)	18.5 (1)	20.5 (3)	27 (4)	29.5 (6)
	Oncomine Precision Assay (cfTNA workflow)	16.5 (1)	17.5 (2)	20 (3)	22.5 (4)
	Oncomine Precision Assay (FFPE workflow)	16.5 (3)	21 (6)	26 (9)	28.5 (12)
	Oncomine TCR Beta-LR Assay	20 (4)	24 (8)	28.5 (12)	31 (16)

# Oncomine Precision Assay on Ion Torrent Genexus System

Maximizes your ability to detect relevant variants

## Curated pan-cancer content



- Mutations, CNVs, and fusion variant types across 50 key genes
- Tumor suppressors, drivers, and resistance variants

## Tissue and plasma samples



- One test, one workflow, multiple sample types
- Maximizes the number of tumors that can be profiled

## Molecular tagging



- Enhanced low-level variant detection
- Key for liquid biopsy testing

## FusionSync™ Detection Technology



- Sensitive and specific—targeted isoform designs
- Novel fusion detection

The content provided herein may relate to products that have not been officially released and is subject to change without notice.

# The Oncomine Precision Assay Content



The Oncomine Precision Assay content is carefully curated to include all relevant targets and targets of emerging importance in precision oncology clinical research..

- 50 genes and 2,769 unique variants
- Mutations (45), CNVs (14), and fusion variants (19),
- Pan cancer span with NSCLC focus
- 218 potential resistance mutations across 22 genes

The content provided herein may relate to products that have not been officially released and is subject to change without notice




# Oncomine Precision Assay Gene Content

DNA hotspots			CNV	Inter-genetic fusions		Intra-genetic fusions
AKT1	ESR1	MAP2K2	ALK	ALK	NTRK2	AR
AKT2	FGFR1	MET	AR	BRAF	NTRK3	BRAF
AKT3	FGFR2	MTOR	CD274	ESR1	NUTM1	EGFR
ALK	FGFR3	NRAS	CDKN2A	FGFR1	RET	MET
AR	FGFR4	NTRK1	EGFR	FGFR2	ROS1	
ARAF	FLT3	NTRK2	ERBB2	FGFR3	RSPO2	
BRAF	GNA11	NTRK3	ERBB3	MET	RSPO3	
CDK4	GNAQ	PDGFRA	FGFR1	NRG1		
CDKN2A	GNAS	PIK3CA	FGFR2	NTRK1		
CHEK2	HRAS	PTEN	FGFR3			
CTNNB1	IDH1	RAF1	KRAS			
EGFR	IDH2	RET	MET			
ERBB2	KIT	ROS1	PIK3CA			
ERBB3	KRAS	SMO	PTEN			
ERBB4	MAP2K1	TP53				

The content provided herein may relate to products that have not been officially released and is subject to change without notice

# Pan-Cancer Clinical Research Application of OPA




The content provided herein may relate to products that have not been officially released and is subject to change without notice

# Oncomine Precision Assay Design Details

The **Oncomine Precision Assay** was designed to detect **2,769** total unique variants (hotspot mutations, CNV genes, and fusion isoforms), with some key statistics:

- **269** TP53 mutations which accounts for approximately 90% of known, recurrent hotspots
- **151** NTRK1/2/3 fusion isoforms across varying partners and breakpoints
- **156** FGFR1/2/3 fusion isoforms across varying partners and breakpoints
- **89** RET fusion isoforms across varying partners and breakpoints
- **3** intra-genetic variants (MET Exon 14 Skip, EGFRviii, and ARv7)
- **218** potential resistance mutations across 22 genes



The content provided herein may relate to products that have not been officially released and is subject to change without notice

# Benefits of Oncomine Fusion Detection Technology



Generally, there are **two** key features for optimal fusion detection:

1. Performance of fusion detection with low input samples / low level transcripts
2. Ability to detect novel fusions for driver genes (e.g. *NTRK* and *FGFR*)


Many similar technologies emphasize #2 above but ignore #1.

With FusionSync™ detection, we address **BOTH #1 and #2**

The content provided herein may relate to products that have not been officially released and is subject to change without notice

# FusionSync Detection Technology

**FusionSync Detection Technology** is a synchronous approach that combines three methods for sensitive, specific, and broad detection of known and novel fusions




\* Available for ALK, FGFR1, FGFR2, FGFR3, NTRK1, NTRK2, NTRK3, and RET fusion drivers



# FusionSync™ Detection Using Genexus and Oncomine Precision Assay

## Detection of ALK and ROS1 fusions


24 NSCLC FFPE samples with orthogonally confirmed ALK or ROS1 fusions were tested:



**100% agreement** in detecting ALK and ROS1 fusions using Oncomine Precision Assay and Genexus

Testing was performed at multiple internal R&D laboratories as part of product verification testing. Samples were run at different plexy levels. Additional verification and validation testing will be performed using the Oncomine Precision Assay and FFPE samples; therefore final performance values can change.

## Detection of RET fusion using **exon tiling imbalance**



Fusion Driver Detected	Predicted Breakpoint Range	Imbalance Score (threshold)	P-Value (threshold)
RET	ex8-ex12	2.12 (> 1.65)	3e-03 (< 0.05)

Testing on LC/2ad cell line that includes CCDC6-RET fusion. Blue line indicates normalized read counts of tiled amplicons across RET, showing a differential expression from the 3' (right) to 5' (left) end of the gene. Collection of grey lines indicate RET expression from fusion wild-type samples, to be used as a baseline comparison to sample RET expression measurement (blue line). Red dotted line indicates predicted break point of RET fusion

# ALK Example: FusionSync™ Detection Using Targeted and Novel Detection Methods

Targeted isoform  
detection



*HIP1-ALK.H21A20*

Isoform read count:  
373

Isoform molecular count:  
3

*KIF5B-ALK.K17A20*

Isoform read count:  
602




Isoform molecular count:  
10

*EML4-ALK.E20A20*

Isoform read count:  
671

Isoform molecular count:  
12

Novel detection  
using imbalance




The content provided herein may relate to products that have not been officially released and is subject to change without notice

# FFPE Samples Successfully Tested Using Oncomine Precision Assay





Uncharacterized FFPE samples were included in internal verification testing using the **Oncomine Precision Assay** and **Genexus Integrated Sequencer**. Results demonstrated successful sequencing of different tumor and biopsy types.



Cancer Type
Bladder
Breast (x2)
Kidney (x2)
Ovarian
Stomach (x2)
Testis
Uterus (x2)

	Cancer Type
FNA	Lung
	Lung
	Lung
CNB	Breast
	Breast
	Lung



*\*NOTE: Samples were not orthogonally tested, therefore variant results are not presented*

**NOTE:** Testing was performed at multiple internal R&D laboratories as part of product verification testing. Samples were run at different plexy levels. Additional verification and validation testing will be performed using the Oncomine Precision Assay and FFPE samples; therefore final performance values can change.

The content provided herein may relate to products that have not been officially released and is subject to change without notice

# Performance of Oncomine Precision Assay on Liquid Biopsy Controls



Following are results from internal verification testing using the **Oncomine Precision Assay** and **Genexus Integrated Sequencer** with liquid biopsy controls.

Sample	Description	Range	Variant Type	# of Variants Per Sample	# of Sample Replicates	Sensitivity	PPV
Internal Liquid Biopsy Control	Synthetic control across 24 genes	0.33% $\pm$ 0.17% SD AF*	SNV	67	32	89.2%	100%
			INDEL	4	32	100%	100%
Internal Liquid Biopsy Control	Synthetic control that contains <i>EGFR</i> amplification	1.16X fold change	CNV	1	32	100%	100%
Internal Liquid Biopsy Control	Synthetic control that contains <i>MET</i> exon 14 skip	1% Tri-Fusion and <i>MET</i> Exon 14 Skip in Total RNA	FUSION	1	32	97.7%	96.4%

**NOTE:** Testing was performed at multiple internal R&D laboratories as part of product verification testing. Samples were run at different plexy levels. Additional verification and validation testing will be performed using the Oncomine Precision Assay and FFPE controls; therefore final performance values can change.  
\* indicates observed allele frequency using Oncomine Precision Assay on Genexus.

The content provided herein may relate to products that have not been officially released and is subject to change without notice

# Key Genes and Variant Types Detected with Cell Lines and FFPE Samples

The following are key variants that were detected using the Oncomine Precision Assay and from cell lines and FFPE samples:

Gene	Variant	Variant Type	Sample Type
<i>EGFR</i>	p.A767_V769dup	Insertion	FFPE
<i>EGFR</i>	p.E746_S752delinsV	Deletion	FFPE
<i>KRAS</i>	G12D	SNV	FFPE
<i>ERBB2 (HER2)</i>	Amplification	CNV	FFPE
<i>ROS1</i>	Fusion	FUSION	FFPE
<i>PTEN</i>	Loss	CNV	Cell Line

**NOTE:** Testing was performed at multiple internal R&D laboratories as part of product verification testing. Samples were run at different plexy levels. Additional verification and validation testing will be performed using the Oncomine Precision Assay and FFPE samples; therefore final performance values can change.

The content provided herein may relate to products that have not been officially released and is subject to change without notice



# Detection of *ALK* and *ROS1* Fusions in FFPE Samples

Sample #	FISH	IHC	OPA on Genexus	Agreement	Imbalance	P-Value	Breakpoint	Concordance
1	ALK	ALK	EML4-ALK	YES	5.313	0.0007	exon15-exon20	YES
2	ALK	ALK	EML4-ALK	YES	3.969	0.0007	exon14-exon15	YES
3	ALK	ALK	EML4-ALK	YES	5.021	0.0007	exon15-exon20	YES
4	ALK	ALK	EML4-ALK	YES	5.313	0.0007	exon15-exon20	YES
5	ALK	ALK	EML4-ALK	YES	5.273	0.0007	exon15-exon20	YES
6	ALK	ALK	EML4-ALK	YES	3.033	0.0007	exon10-exon13	YES
7	ALK	ALK	EML4-ALK	YES	5.189	0.0007	exon15-exon20	YES
8	ALK	ALK	EML4-ALK	YES	5.276	0.0007	exon15-exon20	YES
9	ALK	ALK	EML4-ALK	YES	5.034	0.005	exon15-exon20	YES
10	ALK	ALK	EML4-ALK	YES	5.126	0.0007	exon15-exon20	YES
11	ALK	ALK	EML4-ALK	YES	5.171	0.0007	exon15-exon20	YES
12	ALK	ALK	EML4-ALK	YES	4.757	0.0007	exon15-exon20	YES
13	ALK	ALK	EML4-ALK	YES	4.88	0.001	exon15-exon20	YES
14	ALK	ALK	EML4-ALK	YES	3.688	0.0007	exon15-exon20	YES
15	ALK	ALK	EML4-ALK	YES	5.257	0.0007	exon15-exon20	YES
16	ALK	ALK	HIP1-ALK	YES	5.145	0.0007	exon15-exon20	YES
17	ALK	ALK	EML4-ALK	YES	4.059	0.0007	exon14-exon15	YES
18	ROS1	N/A	CD74-ROS1	YES	N/A	N/A	N/A	N/A
19	ROS1	N/A	EZR-ROS1	YES	N/A	N/A	N/A	N/A
20	ROS1	N/A	EZR-ROS1	YES	N/A	N/A	N/A	N/A
21	ROS1	N/A	EZR-ROS1	YES	N/A	N/A	N/A	N/A
22	ROS1	N/A	CD74-ROS1	YES	N/A	N/A	N/A	N/A
23	ROS1	N/A	CD74-ROS1	YES	N/A	N/A	N/A	N/A
24	ROS1	N/A	SDC4-ROS1	YES	N/A	N/A	N/A	N/A

24 FFPE samples (NSCLC) were tested with the **Oncomine Precision Assay** on **Genexus Integrated Sequencer**.

These samples were previously tested using **FISH** and **IHC** (*ALK* only), and determined to be either *ALK* or *ROS1* fusion positive.

Verification results showed **100% agreement** when using the Oncomine Precision Assay compared to both methods,

In addition, all 17 *ALK* positive samples had concordant detection using the exon tiling imbalance method.

*NOTE: ROS1 does not have exon tiling imbalance capability on Oncomine Precision Assay.*

**NOTE:** Testing was performed at multiple internal R&D laboratories as part of product verification testing. Samples were run at different plexy levels. Additional verification and validation testing will be performed using the Oncomine Precision Assay and FFPE samples; therefore final performance values can change.

# Complete Detection of NTRK Isoforms Using SeraCare FFPE Control

SeraSeq® FFPE NTRK Fusion RNA Reference Material	
Oncomine Precision Assay Detects All NTRK Isoforms	✓ TPM3(7) - NTRK1(10)
	✓ LMNA(11) - NTRK1(11)
	✓ IRF2BP2(1) - NTRK1(10)
	✓ SQSTM1(5) - NTRK1(10)
	✓ TFG(5) - NTRK1(10)
	✓ AFAP1(14) - NTRK2(12)
	✓ NACC2(4) - NTRK2(13)
	✓ QKI(6) - NTRK2(16)
	✓ TRIM24(12) - NTRK2(15)
	✓ PAN3(1) - NTRK2(17)
	✓ ETV6(4) - NTRK3(14)
	✓ ETV6(4) - NTRK3(15)
	✓ ETV6(5) - NTRK3(14)
	✓ ETV6(5) - NTRK3(15)
	✓ BTBD1(4) - NTRK3(14)




Metric	5% Dilution
# of Runs	30
# of NTRK Fusion Isoforms	15
# of Total Datapoints	450
False Negative	4
False Positive	4
<b>Sensitivity</b>	<b>99.1%</b>
<b>PPV</b>	<b>100%</b>

- The SeraCare SeraSeq® FFPE NTRK Fusion RNA Reference Material (01710-1031) contains 15 unique NTRK isoforms (5 NTRK1, 5 NTRK2, and 5 NTRK3).
- A 5% dilution (within a normal background) of the NTRK control was made to reduce the transcript levels of the various fusion isoforms.
- The dilution was tested with multiple replicates using the **Oncomine Precision Assay** and the **Genexus Integrated Sequencer**.
- Verification results showed ability of the assay to detect all 15 NTRK fusions using the targeted isoform designs and measuring the read counts for each isoform.
- Replicate testing demonstrated high sensitivity and specificity at 5% dilutions above 33 fusion reads.

**NOTE:** Testing was performed at multiple internal R&D laboratories as part of product verification testing. Samples were run at different plexy levels, resulting in an average mapped RNA reads of 190,972 for the 10% dilution dataset and an average mapped RNA reads of 691,064 for the 5% dilution dataset. Additional verification and validation testing will be performed using the Oncomine Precision Assay and SeraCare NTRK control; therefore final performance values can change.

# Detection of *ALK*, *FGFR3*, *NTRK1*, and *RET* Fusions Using Exon Tiling Imbalance



Fusion driver	Sample	Predicted breakpoint range	Imbalance score	P-value
<i>ALK</i>	Tri-Fusion Cell Line Mixture with <i>ALK</i> , <i>RET</i> , and <i>ROS1</i> fusions	ex15-ex20	4.74	7e-04
<i>FGFR3</i>	RT4 cell line	ex13-3pUTR	1.74	3.9e-03
<i>NTRK1</i>	KM12 cell line	ex6-ex13	2.58	4.3e-03
<i>RET</i>	Tri-Fusion Cell Line Mixture with <i>ALK</i> , <i>RET</i> , and <i>ROS1</i> fusions	ex8-ex12	2.12	3e-03

The content provided herein may relate to products that have not been officially released and is subject to change without notice

# Performance of Genexus and Oncomine Comprehensive Assay v3

	SNV/INDEL	CNV	Fusions
Sensitivity	98.3%	97.2%	97.6%
PPV	99.0%	N/A	99.6%
Inter-run reproducibility	98.6%	97.8%	95.7%
Intra-run reproducibility	98.6%	96.7%	95.9%

*SNV and INDEL performance was assessed using Acrometrix Hotspot Control and SeraCare Seraseq Tri-Level Mutation DNA Mix*

*CNV performance was assessed using SeraCare CNV controls*

*Fusion performance was assessed using Horizon HD789 and SeraCare Fusion Mix v3, Horizon HD784 ALK, and LC/2 ad cell line*

The content provided herein may relate to products that have not been officially released and is subject to change without notice

# A New world of NGS Means a New Day for your Lab

Now is the time to start with NGS with...



- ✓ One day Ion Torrent Genexus workflow\*
- ✓ Complete automation from specimen to report and only 10 min hands on time\*
- ✓ Oncomine Precision Assay
- ✓ No need for batching
- ✓ First in class implementation support

..and get key 50 gene NGS profile result in the same time as your IHC for comprehensive profile

\*Specimen-to-report workflow will be available after the Ion Torrent™ Genexus™ Purification System and integrated reporting capabilities are added in 2020.  
The content provided herein may relate to products that have not been fully validated by Thermo Fisher Scientific and is subject to change without notice.



*The **Ion Torrent™ Genexus™** is currently Research Use Only, however “Thermo Fisher Scientific intent to seek regulatory marketing authorization of the system so that it can be potentially made available in every clinical setting. Additionally we also plan to develop and seek approval for a broad portfolio of diagnostic assays in oncology. The **Ion Torrent™ Genexus™** has been developed with the intent to shift the cancer testing paradigm in the future, that is what we worked for and will keep on working for”.*



# Thank you Q&A

**ThermoFisher**  
SCIENTIFIC

The world leader in serving science



**For Research Use Only** © 2019 Thermo Fisher Scientific Inc. All rights reserved.  
All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Approval: 5748

The world leader in serving science





# Kit Configuration, Throughput and Pricing

**ThermoFisher**  
SCIENTIFIC

The world leader in serving science

# Oncomine Precision Assay Throughput on Genexus



Assay		1 Lane	2 Lanes	3 Lanes	4 Lanes
Max Number of Samples	Oncomine Precision Assay (for DNA & RNA workflow)	4	8	12	16
	Oncomine Precision Assay (for cfTNA workflow)	1	2	3	4



# OPA Content – Progression on Other Oncomine Assays



Swap genes, emphasis on potential for future targeted therapies research

**Examples:**

*NRG1* fusions  
*PTEN* loss  
*AKT2/3* mutations

Increased fusion isoforms and exon tiling imbalance

*NTRK* isoforms **31** (OFA)→**150** (OPA)


**Imbalance:** *ALK*, *FGFR1*, *FGFR2*, *FGFR3*,  
*NTRK1*, *NTRK2* *NTRK3*, *RET*

Resistance mutations

*ALK*, *AR*, *EGFR*, *ERBB2*, *ERBB3*, *ESR1*,  
*FGFR1*, *FGFR2*, *FLT3*, *IDH1*, *IDH2*, *KIT*,  
*MAP2K1*, *MAP2K2*, *MET*, *NTRK1*, *NTRK2*,  
*NTRK3*, *PDGFRA*, *RET*, *ROS1*, *SMO*

The content provided herein may relate to products that have not been officially released and is subject to change without notice

# Assay Roadmap On Genexus





# A New world of NGS Means a New Day for your Lab

Now is the time to start with NGS with...



- ✓ One day Ion Torrent Genexus workflow\*
- ✓ Complete automation from specimen to report and only 10 min hands on time\*
- ✓ Oncomine Precision Assay
- ✓ No need for batching
- ✓ First in class implementation support

..and get key 50 gene NGS profile result in the same time as your IHC for comprehensive profile

\*Specimen-to-report workflow will be available after the Ion Torrent™ Genexus™ Purification System and integrated reporting capabilities are added in 2020.  
The content provided herein may relate to products that have not been fully validated by Thermo Fisher Scientific and is subject to change without notice.

# First Experience from the “New World”

**The Pathologist - Webinar**

**5<sup>th</sup> December 2019**

*José Luis Costa*

([jcosta@ipatimup.pt](mailto:jcosta@ipatimup.pt))





# Ipatimup – Porto - Portugal



- Biggest research in health institute in Portugal (1250 researchers);
- Cancer, Neurosciences and Host-Pathogen interactions research lines




- Founding member of i3S
- Leading cancer research institute in Portugal
- Founding partner of **Porto Comprehensive Cancer Center**






# NGS clinical research timeline....



## Genexus System



# First Experience from IPATIMUP, Portugal

## Commercial control samples

- Horizon and SeraCare FFPE and cfDNA reference material

## Clinical research samples

- Lung cancer tissue and liquid biopsy research samples
- Previously characterize on an Ion S5XL system


# First Experience from IPATIMUP, Portugal

Genexus Software | Ion Torrent

Samples ▾ Runs ▾ Monitor ▾ Results ▾   




Results > Run Results > View Results

OPA - GX5 - Liquid Biopsy TN... Assay: OPA - GX5 - Liquid Biopsy TNA - EA2-20190827  
Run Plan : Run02\_27AUG19\_OPA\_IPEAPh2



# First Experience from IPATIMUP, Portugal


Genexus Software | Ion Torrent

Samples ▾ Runs ▾ Monitor ▾ Results ▾   

Results > Run Results > View Results

OPA - GX5 - Liquid Biopsy TN... Assay: OPA - GX5 - Liquid Biopsy TNA - EA2-20190827  
Run Plan : Run02\_27AUG19\_OPA\_IPEAPh2

Loading Density



89.21%  
Loading

5.84G  
Total Bases

58,534,613  
Final Reads

99.06%  
Mean Raw Accuracy 1x

Loading 89.21%

120,856,893

Wells with ISPs / Total Addressable Wells  
120,856,893 / 135,481,643

Empty Wells 10.79%

Enrichment 99.29%

119,996,579

Wells with Live ISPs / Wells with ISPs  
119,996,579 / 120,856,893

Library 99.99%

119,985,204

Wells with Library ISPs / Wells with live ISPs  
119,985,204 / 119,996,579

Wells with Control ISPs 0.01%  
11,375 / 119,996,579

Final Reads 48.78%

58,534,613

Final Reads / Wells with Library ISPs  
58,534,613 / 119,985,204

■ Polyclonal 32.07%  
38,475,706 / 119,985,204  
■ Low Quality 18.25%  
21,901,337 / 119,985,204  
■ Adapter Dimer 0.02%  
29,212 / 119,985,204

Run Samples ?

Sample Name	Nucleic Acid Type	Barcode	Total Reads	Mean Read Length	>= Q20 Bases	Uniformity	Read Length Histogram
cFTNA_1_27AUG19_IPEAPh2	TNA	IonHDdual_0125	13,150,825	101 bp	1,150,409,567	Not Available	
cFTNA_2_27AUG19_IPEAPh2	TNA	IonHDdual_0126	13,700,575	100 bp	1,192,236,490	Not Available	
cFTNA_3_27AUG19_IPEAPh2	TNA	IonHDdual_0127	16,994,843	96 bp	1,420,277,918	Not Available	


# First Experience from IPATIMUP, Portugal

Genexus Software | Ion Torrent

Samples ▾ Runs ▾ Monitor ▾ Results ▾   

Results > Run Results > View Results

OPA - GX5 - Liquid Biopsy TN... Assay: OPA - GX5 - Liquid Biopsy TNA - EA2-20190827  
Run Plan : Run02\_27AUG19\_OPA\_IPEAPh2






# First Experience from IPATIMUP, Portugal

Genexus Software | Ion Torrent

Samples ▾ Runs ▾ Monitor ▾ Results ▾   

Results > Run Results > View Results




OPA - GX5 - Liquid Biopsy TN... Assay: OPA - GX5 - Liquid Biopsy TNA - EA2-20190827  
Run Plan : Run02\_27AUG19\_OPA\_IPEAPh2






# First Experience from IPATIMUP, Portugal

Genexus Software | Ion Torrent

Samples ▾ Runs ▾ Monitor ▾ Results ▾   

Results > Run Results > View Results

<< < cftNA\_1\_27AUG19\_IPEAPh2 > >> QC Status:  Assay: OPA - GX5 - Liquid Biopsy TNA - EA2-20190827  
Run Plan: Run02\_27AUG19\_OPA\_IPEAPh2

## Sample Details

Sample Name:	cftNA_1_27AUG19_IPEAPh2	Collection Date:	23 AUG 2019
Gender:	Unknown	Sample Type:	cftNA
Disease Category:	Cancer	Cancer Type:	Unknown Primary Origin
Cancer Stage:	Unknown	% Cellularity:	null

## Metrics

Average Base Coverage Depth:	Not Available
Uniformity Of Base Coverage:	Not Available
% Base Reads On Target:	Not Available
Median Molecular Coverage:	2253
Median Read Coverage:	26148

## Variant Summary

A default filter has been applied. Go to [SNVs/Indels](#), [Fusions](#), [CNVs](#) pages to remove or modify variant filter.

Filter Chain Applied: [Variant Matrix tab Summary](#)

### SNVs/Indels

6 Detected

Gene	AA Change	Mol Freq %	Oncomine Variant Class
ERBB3	p.E332K	3.7131	Hotspot
KRAS	p.A59T	0.1335	Hotspot
KRAS	p.G12A	2.6452	Hotspot
MET	p.?	0.263	Hotspot
TP53	p.R248L	1.4074	Hotspot
TP53	p.S241C	0.1141	Hotspot

### Fusions

5 Detected

Oncomine Driver Gene	Evidence Level
MET	Targeted Isoforms
ALK	Targeted Isoforms
BRAF	Targeted Isoforms
RET	Targeted Isoforms
ROS1	Targeted Isoforms




### CNVs

2 Detected


Gene	Gain/Loss	Oncomine Variant Class
EGFR	↑	Amplification
AR	↓	

# First Experience from IPATIMUP, Portugal

Genexus Software | Ion Torrent

Samples ▾ Runs ▾ Monitor ▾ Results ▾   

Results > Run Results > View Results

<< < cftNA\_1\_27AUG19\_IPEAPh2 > >> QC Status:  Assay: OPA - GX5 - Liquid Biopsy TNA - EA2-20190827  
Run Plan: Run02\_27AUG19\_OPA\_IPEAPh2

## Sample Details

Sample Name:	cftNA_1_27AUG19_IPEAPh2	Collection Date:	23 AUG 2019
Gender:	Unknown	Sample Type:	cftNA
Disease Category:	Cancer	Cancer Type:	Unknown Primary Origin
Cancer Stage:	Unknown	% Cellularity:	null

## Metrics

Average Base Coverage Depth:	Not Available
Uniformity Of Base Coverage:	Not Available
% Base Reads On Target:	Not Available
Median Molecular Coverage:	2253
Median Read Coverage:	26148

## Variant Summary

A default filter has been applied. Go to [SNVs/Indels](#), [Fusions](#), [CNVs](#) pages to remove or modify variant filter.

Filter Chain Applied: Variant Matrix tab Summary

### SNVs/Indels

6 Detected

Gene	AA Change	Mol Freq %	Oncomine Variant Class
ERBB3	p.E332K	3.7131	Hotspot
KRAS	p.A59T	0.1335	Hotspot
KRAS	p.G12A	2.6452	Hotspot
MET	p.?	0.263	Hotspot
TP53	p.R248L	1.4074	Hotspot
TP53	p.S241C	0.1141	Hotspot

### Fusions

5 Detected

Oncomine Driver Gene	Evidence Level
MET	Targeted Isoforms
ALK	Targeted Isoforms
BRAF	Targeted Isoforms
RET	Targeted Isoforms
ROS1	Targeted Isoforms

### CNVs

2 Detected

Gene	Gain/Loss	Oncomine Variant Class
EGFR	↑	Amplification
AR	↓	

# First Experience from IPATIMUP, Portugal

## SNVs/InDels

Sample ID	Gene	Variant	Molecular frequency (%)	
			Run 1	Run 2
FFPE1	AKT1	p.?	4,2	4,9
FFPE1	BRAF	p.V600E	15,5	19,3
FFPE1	CDKN2A	p.R58*	nd	9,3
FFPE1	CTNNB1	p.S33Y	5,9	6,2
FFPE1	CTNNB1	p.S45del	4,8	4,3
FFPE1	EGFR	p.A767_V769dup	3,5	5,1
FFPE1	EGFR	p.E746_A750del	3,1	nd
FFPE1	EGFR	p.G719S	7,4	7,1
FFPE1	GNA11	p.Q209L	3,9	2,7
FFPE1	KRAS	p.G13D	2,8	3,7
FFPE1	MAP2K1	p.Q56P	3,8	3,4
FFPE1	PIK3CA	p.E545K	3,8	3,5
FFPE1	PIK3CA	p.H1047R	19,8	17,6
FFPE1	TP53	p.S241C	9,0	5,8
FFPE1	TP53	p.S241F	9,0	10,2
FFPE2	AKT3	p.Q78K	nd	2,2
FFPE2	BRAF	p.V600E	64,3	69,7
FFPE2	EGFR	p.E746_A750del	3,5	3,6
FFPE2	EGFR	p.G719S	5,4	6,3
FFPE2	EGFR	p.L858R	3,9	5,1
FFPE2	EGFR	p.L861Q	3,9	3,4
FFPE2	EGFR	p.T790M	4,9	4,1
FFPE2	PIK3CA	p.H1047R	50,8	54,6
FFPE3	CTNNB1	p.S33Y	51,6	51,0
FFPE3	EGFR	p.G719S	31,5	32,6
FFPE3	KRAS	p.A146T	3,8	3,7
FFPE3	KRAS	p.G12D	5,9	4,4
FFPE3	KRAS	p.G13D	4,1	4,0
FFPE3	KRAS	p.Q61H	3,9	3,1
FFPE3	MAP2K1	p.Q56P	49,4	52,2
FFPE3	NRAS	p.G12V	4,0	6,0
FFPE3	NRAS	p.Q61K	3,6	3,8
FFPE5	BRAF	p.V600E	15,3	14,3
FFPE5	EGFR	p.E746_A750del	13,1	10,8
FFPE5	EGFR	p.L858R	14,6	16,1
FFPE5	ERBB2	p.Y772_A775dup	15,9	14,8
FFPE5	IDH1	p.R132C	14,3	13,3
FFPE5	KRAS	p.G12D	14,2	14,5
FFPE5	PIK3CA	p.E542K	16,3	16,4

## Fusions

Sample ID	Driver Gene	Run 1	Run 2
FFPE2	ALK	Detected	Detected
FFPE2	BRAF	Detected	Detected
FFPE2	NTRK1	Detected	Detected
FFPE2	NTRK2	Detected	Detected
FFPE2	NTRK3	Detected	Detected
FFPE5	FGFR3	Detected	Detected

## CNVs

Sample ID	Gene	Run 1	Run 2
FFPE1	CDKN2A	loss	loss
FFPE1	MET	gain	gain
FFPE5	MET	gain	gain
FFPE5	AR	loss	Loss
FFPE5	FGFR3	loss	nd

- Different variant types can be detected;

# First Experience from IPATIMUP, Portugal

## SNVs/InDels

Sample ID	Gene	Variant	Molecular frequency (%)	
			Run 1	Run 2
FFPE1	AKT1	p.?	4,2	4,9
FFPE1	BRAF	p.V600E	15,5	19,3
FFPE1	CDKN2A	p.R58*	nd	9,3
FFPE1	CTNNB1	p.S33Y	5,9	6,2
FFPE1	CTNNB1	p.S45del	4,8	4,3
FFPE1	EGFR	p.A767_V769dup	3,5	5,1
FFPE1	EGFR	p.E746_A750del	3,1	nd
FFPE1	EGFR	p.G719S	7,4	7,1
FFPE1	GNA11	p.Q209	3,9	2,7
FFPE1	KRAS	p.G13D	2,8	3,7
FFPE1	MAP2K1	p.Q56P	3,8	3,4
FFPE1	PIK3CA	p.E542K	3,8	3,5
FFPE1	PIK3CA	p.H1047R	19,8	17,6
FFPE1	TP53	p.S241C	9,0	5,8
FFPE1	TP53	p.S241F	9,0	10,2
FFPE2	AKT3	p.Q78K	nd	2,2
FFPE2	BRAF	p.V600E	64,3	69,7
FFPE2	EGFR	p.E746_A750del	3,5	3,6
FFPE2	EGFR	p.G719S	5,4	6,3
FFPE2	EGFR	p.L858R	3,9	5,1
FFPE2	EGFR	p.L861Q	3,9	3,4
FFPE2	EGFR	p.T792M	4,9	4,1
FFPE2	PIK3CA	p.H1047R	50,8	54,6
FFPE3	CTNNB1	p.S33Y	51,6	51,0
FFPE3	EGFR	p.G719S	31,5	32,6
FFPE3	KRAS	p.A146T	3,8	3,7
FFPE3	KRAS	p.G12D	5,9	4,4
FFPE3	KRAS	p.G13D	4,1	4,0
FFPE3	KRAS	p.Q61H	3,9	3,1
FFPE3	MAP2K1	p.Q56P	49,4	52,2
FFPE3	NRAS	p.G12V	4,0	6,0
FFPE3	NRAS	p.Q61K	3,6	3,8
FFPE5	BRAF	p.V600E	15,3	14,3
FFPE5	EGFR	p.E746_A750del	13,1	10,8
FFPE5	EGFR	p.L858R	14,6	16,1
FFPE5	ERBB2	p.Y772_A775dup	15,9	14,8
FFPE5	IDH1	p.R132C	14,3	13,3
FFPE5	KRAS	p.G12D	14,2	14,5
FFPE5	PIK3CA	p.E542K	16,3	16,4

## Fusions

Sample ID	Driver Gene	Run 1	Run 2
FFPE2	ALK	Detected	Detected
FFPE2	BRAF	Detected	Detected
FFPE2	NTRK1	Detected	Detected
FFPE2	NTRK2	Detected	Detected
FFPE2	NTRK3	Detected	Detected
FFPE5	FGFR3	Detected	Detected

## CNVs

Sample ID	Gene	Run 1	Run 2
FFPE1	CDKN2A	loss	loss
FFPE1	MET	gain	gain
FFPE5	MET	gain	gain
FFPE5	AR	loss	Loss
FFPE5	FGFR3	loss	nd

- Different variant types can be detected;
- Reproducibility high between runs;
- Curiously, deviation of molecular frequency is bigger for higher molecular frequencies (MF<10% SD:0.6 – MF>10% SD:1.4);

# First Experience from IPATIMUP, Portugal

## SNVs/InDels

Sample ID	Gene	Variant	Molecular frequency (%)	
			Run 1	Run 2
FFPE1	AKT1	p.?	4,2	4,9
FFPE1	BRAF	p.V600E	15,5	19,3
FFPE1	CDKN2A	p.R58*	nd	9,3
FFPE1	CTNNB1	p.S33Y	5,9	6,2
FFPE1	CTNNB1	p.S45del	4,8	4,3
FFPE1	EGFR	p.A767_V769dup	3,5	5,1
FFPE1	EGFR	p.E746_A750del	3,1	nd
FFPE1	EGFR	p.G719S	7,4	7,1
FFPE1	GNA11	p.Q209L	3,9	2,7
FFPE1	KRAS	p.G13D	2,8	3,7
FFPE1	MAP2K1	p.Q56P	3,8	3,4
FFPE1	PIK3CA	p.E545K	3,8	3,5
FFPE1	PIK3CA	p.H1047R	19,8	17,6
FFPE1	TP53	p.S241C	9,0	5,8
FFPE1	TP53	p.S241F	9,0	10,2
FFPE2	AKT3	p.Q78K	nd	2,2
FFPE2	BRAF	p.V600E	64,3	69,7
FFPE2	EGFR	p.E746_A750del	3,5	3,6
FFPE2	EGFR	p.G719S	5,4	6,3
FFPE2	EGFR	p.L858R	3,9	5,1
FFPE2	EGFR	p.L861Q	3,9	3,4
FFPE2	EGFR	p.T790M	4,9	4,1
FFPE2	PIK3CA	p.H1047R	50,8	54,6
FFPE3	CTNNB1	p.S33Y	51,6	51,0
FFPE3	EGFR	p.G719S	31,5	32,6
FFPE3	KRAS	p.A146T	3,8	3,7
FFPE3	KRAS	p.G12D	5,9	4,4
FFPE3	KRAS	p.G13D	4,1	4,0
FFPE3	KRAS	p.Q61H	3,9	3,1
FFPE3	MAP2K1	p.Q56P	49,4	52,2
FFPE3	NRAS	p.G12V	4,0	6,0
FFPE3	NRAS	p.Q61K	3,6	3,8
FFPE5	BRAF	p.V600E	15,3	14,3
FFPE5	EGFR	p.E746_A750del	13,1	10,8
FFPE5	EGFR	p.L858R	14,6	16,1
FFPE5	ERBB2	p.Y772_A775dup	15,9	14,8
FFPE5	IDH1	p.R132C	14,3	13,3
FFPE5	KRAS	p.G12D	14,2	14,5
FFPE5	PIK3CA	p.E542K	16,3	16,4

## Fusions

Sample ID	Driver Gene	Run 1	Run 2
FFPE2	ALK	Detected	Detected
FFPE2	BRAF	Detected	Detected
FFPE2	NTRK1	Detected	Detected
FFPE2	NTRK2	Detected	Detected
FFPE2	NTRK3	Detected	Detected
FFPE5	FGFR3	Detected	Detected

## CNVs

Sample ID	Gene	Run 1	Run 2
FFPE1	CDKN2A	loss	loss
FFPE1	MET	gain	gain
FFPE5	MET	gain	gain
FFPE5	AR	loss	Loss
FFPE5	FGFR3	loss	nd

- Different variant types can be detected;
- Reproducibility high between runs;
- Curiously, deviation of molecular frequency is bigger for higher molecular frequencies (MF<10% SD:0.6 – MF>10% SD:1.4);

Accuracy of 98% (95%CI 94% to 99%);

Sensitivity of 96 % (95%CI 89% to 99%)

Specificity of 100% (95%CI 94% to 100%).

# First Experience from IPATIMUP, Portugal

## Commercial control samples

- Horizon and SeraCare FFPE and cfDNA reference material

## Clinical research samples

- Lung cancer tissue and liquid biopsy research samples
- Previously characterize on an Ion S5XL system



# Tissue biopsies – RNA and DNA

Lung cancer FFPE sample were sequenced using Colon and Lung or Lung Fusion panel on Ion S5XL system and the Oncomine Precision Assay on Genexus instrument

Variants covered by both panels were detected in both systems at similar allelic frequencies

Additional variants were detected using the Oncomine Precision Assay

Sample 1 presented additional p.T790M mutation that had not been previously identified (intra-sample heterogeneity?)



FFPE	Gene	AA Change	Genexus	Ion S5XL
Sample 1	<b>EGFR</b>	<b>p.L858R</b>	<b>43,8</b>	<b>31,3</b>
	TP53	p.R248Q	1,8	nd
	<b>EGFR</b>	<b>p.T790M</b>	<b>1,6</b>	<b>nd</b>
	ALK	p.A1200V	1,5	nd
Sample 2	<b>EGFR</b>	<b>p.L747_P753delinsS</b>	<b>28,6</b>	<b>20,1</b>
	TP53	p.C176R	3,9	nd
	RET	p.G810S	1,6	nd
	EGFR	p.P848L	1,6	nd
	FGFR3	p.R399C	0,3	nd
Sample 3	<b>ALK</b>	<b>fusion</b>	<b>detected</b>	<b>detected</b>
Sample 4	<b>BRAF</b>	<b>p.V600E</b>	<b>47,8</b>	<b>36,3</b>
	TP53	p.R175C	13,8	nd
	FGFR2	p.A648T	6,9	nd
	EGFR	p.R836C	6,2	nd
	ERBB3	p.V104M	4,4	nd
	PDGFRA	p.T849C	4,2	nd
	PIK3CA	p.V344M	3,7	nd
	MAP2K1	p.K57N	2,8	nd
	GNAQ	p.R183Q	2,1	nd
	EGFR	p.V769M	1,8	nd
Sample 5	<b>KRAS</b>	<b>p.G12D</b>	<b>58,3</b>	<b>44,7</b>
	TP53	p.C176Y	31,4	nd
	FGFR3	p.R399C	3,6	nd
	CDKN2A	p.R58Q	3,1	nd
	RET	p.R912Q	2,5	nd
	BRAF	p.D594N	2,3	nd

nd – not detected

# Liquid biopsies – Lung cancer plasma sample at recurrence

9AM  
10min  
hands on

Plasma  
collection


## 1. Plasma collection

BD Vacutainer PPT  
(K2EDTA)



## 2. cfTNA isolation

MagMax cfTNA  
extraction kit



11AM  
90min  
hands on

cfTNA  
isolation

## 3. NGS variant identification



Genexus

+

Oncomine  
Precision Assay

Gene	AA Change	Mol. Freq. %
EGFR	p.E746_A750del	1,59
EGFR	p.T790M	1,78

4PM  
10min  
hands on

Variant  
detection

## 4. OKR Report

ipatimup

IPATIMUP  
Rua Júlio Amaral de Carvalho 45  
4200-135 Porto  
Portugal

Patient ID: Tracking No: Date: 31 May 2019 1 of 2

Name: Primary Tumor Site:  
Date of Birth: Sample Type:  
Gender: Sample ID:  
Case Number: MD/Hospital:  
Entrance Date:

Sample Cancer Type: Melanoma

Report Highlights  
1 Clinically Significant Biomarkers  
7 Therapies Available  
27 Clinical Trials

Biomarker Result

Relevant Therapies Relevant Therapies

- Variant Result
- Interpretation

9AM  
10min  
hands on

Report

Sample

< 24h

Report

# Genexus implementation scenarios

## Stand-alone

All-in-one



Genexus

High volume  
Routine assays

Occasional assays  
Short TAT assays

# Genexus implementation scenarios

## Stand-alone

All-in-one



Genexus

High volume  
Routine assays

Occasional assays  
Short TAT assays

## Bundle our strategy

Workhorse



Ion Chef



Ion S5XL

High volume  
Routine assays

Fast track



Genexus

Occasional assays  
Short TAT assays







Joana Reis  
Sandra Coelho  
José Carlos Machado



Venceslau Hespanhol  
Gabriela Fernandes  
Fatima Carneiro  
Conceição Souto Moura



Sohaib Qureshi  
Jussi Vanhatalo  
Ian Grinsell  
Rosella Petraroli  
Andy Felton



## Funding



Thermo Fisher Scientific and its affiliates are not endorsing, recommending, or promoting any use or application of Thermo Fisher Scientific products presented by third parties during this seminar. Information and materials presented or provided by third parties are provided as-is and without warranty of any kind, including regarding intellectual property rights and reported results. Parties presenting images, text and material represent they have the rights to do so.

Speaker was provided honorarium by Thermo Fisher Scientific for this presentation