

# Comprehensive and simultaneous genomic profiling

The OncoPrint™ Comprehensive Assay Plus for broad pan-cancer profiling across multiple biomarkers and cancer tissue types

Designed for the detection of known and novel biomarkers associated with targeted and immuno-oncology therapies, this new solution for the testing of solid tumor samples covers over 500 genes associated with a variety of tumor types, offering myriad potential pan-cancer research implications.

The Ion Torrent™ OncoPrint™ Comprehensive Assay Plus enables simultaneous analysis of both DNA and RNA in a single workflow, and fits seamlessly into existing workflows for Ion Chef™ and Ion GeneStudio™ S5 systems, providing relevant variant data from integrated informatics.

## Key features of the OncoPrint™ Comprehensive Assay Plus include:

- Over 500 unique genes, including driver genes across solid tumors
- Comprehensive and simultaneous genomic profiling (from DNA and RNA) for research of key targeted and immunotherapy biomarkers
- All classes of variants: SNVs, indels, **CNVs**, fusions, and splice variants
- Immuno-oncology content: microsatellite instability (MSI) and tumor mutational burden (TMB)
- **Sensitive** targeted fusion **detection**
- **Low sample input** (tested with sample inputs down to 1 ng of DNA/RNA or a single slide)
- **Very low QNS**, enabling up to 99% sequencing success [1, 2]
- Fully integrated analysis and reporting
- Fast turnaround time, driven by an automated workflow that enables you to go from sample to report in days

The OncoPrint Comprehensive Assay Plus is part of OncoPrint™ Solutions, a growing menu of assays designed for the way you work. With options ranging from a small, specialized two-gene assay for *BRCA* to a large >500 gene assay for multibiomarker profiling across multiple cancer types, OncoPrint Solutions can help you choose the best solution for your specific oncology research needs.

## Streamlined NGS workflow

### OncoPrint informatics

Our streamlined Ion Torrent™ OncoPrint™ informatics workflow is an integrated, internally developed sample-to-report solution for data analysis that takes you from initial sequence analysis of many variants through annotation of relevant biomarkers associated with cancer drivers to a final report that is fully customizable and can be translated into multiple languages. A sample report from the OncoPrint Comprehensive Assay Plus workflow is at right.

### References

1. Anna-Lena Volckmar, Jonas Leichsenring, Martina Kirchner et al. (2019) Combined targeted DNA and RNA sequencing of advanced NSCLC in routine molecular diagnostics: Analysis of the first 3,000 Heidelberg cases. *Int J Cancer*. <https://doi.org/10.1002/ijc.32133>.
2. Chih-Jian Lih, Robin D. Harrington, David J. Sims et al. (2017) Analytical validation of the next-generation sequencing assay for a nationwide signal-finding clinical trial. *J Mol Diagn* 19(2):313-327. <https://doi.org/10.1016/j.jmoldx.2016.10.007>.

Find out more at [thermofisher.com/ocaplus](https://thermofisher.com/ocaplus)

## OncoPrint Comprehensive Assay Plus: Assay performance

**A** Assay performance on reference control (Thermo Scientific™ AcroMetrix™ Oncology Hotspot Control)


Variant type	TP	FP	FN	Sensitivity	PPV
SNV	2,181	4	33	98.51%	99.82%
Indel	62	2	1	98.41%	96.88%


**B** Assay performance on CNV

Variant type	TP	FN	TN	Sensitivity	Specificity
CNV gain	34	1	12	97.1%	100%
CNV loss	17	2	12	89.4%	

**Figure 1. Distributions of total and amplicon-mapped DNA reads (A) show excellent on-target rates across 243 FFPE samples from a diverse set of tissue types.** A high proportion of amplicon-assigned reads reflect high assay specificity in target amplification. A median of ~35 million reads per sample provide >2,400 average read depth. Reads were distributed across amplicons with >95% uniformity and at least 500x base coverage for >95% of the targets (B) to ensure detection of variants with low allele frequency (>5%).



## OncoPrint Comprehensive Assay Plus: TMB performance



**Figure 2. TMB performance of OncoPrint Comprehensive Assay Plus.** Whole-exome sequencing (WES) is viewed as the gold standard for TMB quantitation. *In silico* analysis against WES was performed to characterize the TMB performance of the OncoPrint Comprehensive Assay Plus. (A) A scatter plot shows high correlation between the targeted assay (y-axis) and WES (x-axis) mutation counts, which was downloaded from TCGA MC3<sup>3</sup>. (B) A second comparison was conducted against a commercially available assay with real-world samples, and also shows high correlation in TMB determination.