

Characterization Of The Enhanced Enzyme For Isothermal Amplification – Invitrogen™ Lyo-ready Bst DNA Polymerase

J Zujienė, U Krukonytė, A Klimentov, L Gabartienė, V Dubovskaja, A Alminaitė, R Sukackaitė
Thermo Fisher Scientific Baltics, Vilnius, Lithuania

Abstract

Note: the abstract has been revised since submission to expand LR Bst active temperatures and exclude Xylan for readability.

Purpose: Characterization of the Invitrogen™ Lyo-ready Bst DNA Polymerase, in-vitro enhanced enzyme, Large Fragment of the Bst DNA polymerase (LR Bst).

Methods: Polymerase activity was tested using radiolabelled dNTPs at 60°C. This temperature was varied during active temperature range testing. Functional testing was carried out using loop-mediated isothermal amplification (LAMP). The method was also used to test polymerase inhibitor resistance by spiking the reaction with inhibitors.

Results: Activity testing revealed that LR Bst active temperature range is 15-79°C. At 37°C, LR Bst has retained ~45% of its activity compared to the activity at 60°C. The increase in activity was linear with regards to increasing temperature. LAMP testing also showed that the Invitrogen Lyo-ready Bst DNA Polymerase has high salt tolerance. The addition of 150 mM KCl should not impact reaction time significantly and, whereas the addition of 125 mM KCl to most LAMP reactions using competitor Bst polymerases inhibited DNA amplification, our LR Bst was able to withstand up to 200 mM additional KCl to the optimized buffer. Inhibitor tolerance was tested further by adding various common amplification inhibitors. Adding up to 5% urea to the reaction failed to significantly increase time to reaction when using the LR Bst. In addition, the Invitrogen Lyo-ready Bst DNA Polymerase showed consistently faster amplification speed with titres of other inhibitors, such as humic acid, or ethanol, than competitor Bst polymerases. Overall, the Invitrogen Lyo-ready Bst DNA Polymerase's ammonium salt and glycerol-free composition, wide active temperature range, and high speed brings an extreme advantage to diagnostic assay developers designing assays for field or clinical applications.

Introduction

Most pathogens that infect humans are of animal origin¹. While a number of new zoonotic diseases have been recorded in the past century, the ongoing antibiotic resistance continues to pose the risk of novel zoonoses emerging². The COVID-19 pandemic has shown that, when it comes to rapidly spreading novel diseases, preventative measures through quick and reliable diagnostics are key if therapeutic options are unavailable³.

The current gold standard to carry out nucleic acid amplification in molecular diagnostics is the polymerase chain reaction (PCR)⁴. Though it shows great specificity and sensitivity, the time and resources needed for the reaction may not be sustainable when a rapid response is required. During the COVID-19 pandemic, where the first-line defense was largely based on diagnosis, the world has seen a significant increase to attempt adapting isothermal nucleic acid amplification techniques (NAATs) to molecular diagnostics⁵. In fact, multiple SARS-CoV-2 LAMP tests have been approved by the U.S. Food & Drug Administration, some of which may be used at home⁶. Isothermal nucleic acid amplification allows for sample detection at constant temperatures, which can be applicable in a wider range of settings than PCR. In addition, whole assays can be lyophilized to provide room temperature stability and easy use at reconstruction⁷. Bst polymerases are some of the most versatile polymerases for LAMP; their strand-displacement feature combined with thermostability allows for application in a range of NAATs. Thus, to address demands of molecular diagnostic assay developers using NAATs, Thermo Fisher Scientific has developed an evolved LR Bst the Invitrogen Lyo-ready Bst DNA Polymerase.

Materials and methods

Polymerase activity testing

Polymerase activity testing for the Invitrogen Lyo-ready Bst DNA Polymerase (Thermo Fisher Scientific) was conducted using an in-house test with radiolabelled nucleotides. One unit is the amount of enzyme that will incorporate 10 pmol of dNTP into a polynucleotide fraction at 65°C in 30 minutes. Incubation temperatures were varied for temperature activity testing.

LAMP


All LAMP and RT-LAMP reactions were assembled according to vendor recommendations (Thermo Fisher Scientific for Invitrogen Lyo-ready Bst DNA Polymerase, New England Biolabs for Bst2.0 WarmStart DN A Polymerase and Bst 3.0 DNA Polymerase, Oligene for GspSD2.0 3 U of Bst, LF were used per (RT-)LAMP reaction. Reaction inhibitors were mixed in-house. They were mixed with sample before addition to the reaction. Targets: Adenovirus1 (ADV1) DNA (Viroset), Mycoplasma pneumoniae PH strain of Eaton Agent (NCTC 10119) (ATCC); Synthetic SARS-CoV-2 RNA control 1, Measles virus RNA control (Twist Bioscience). Primers were ordered from Metabion. Adv141-specific hexon targeting primers designed based on Ziroz et al. (2015)⁸, M. pneumoniae primers from Arfaatabar et al. (2019)⁹. Measles primers designed based on Fujino et al. (2005)¹⁰. SARS-CoV-2 primers targeting designed based on Alekseenko et al. (2021)¹¹.

Results

Invitrogen Lyo-ready Bst DNA Polymerase activity testing

Bst DNA Polymerases may be used in a number of isothermal amplification reactions at various reaction temperatures, based on primers and the thermostability of other enzymes. Thus, Invitrogen Lyo-ready Bst DNA Polymerase activity was first measured at temperatures 35-75°C to reflect its applicability in whole genome amplification (WGA), rolling circle amplification (RCA), and LAMP^{12, 13}. Then, down to 15°C due to promising activity results and possible applicability in low temperature reactions. This was achieved by measuring the LR Bst activity at different temperatures and normalizing the results to that of activity at 60°C (Figure 1). Overall, we observed a stepwise increase in activity from 60 to 65°C (22%) suggesting an optimal temperature range for the LR Bst of 65 to 70°C (20%). A further 11% increase in activity from 70 to 75°C was also observed, which could indicate entering suboptimal conditions for LR Bst. Activity at even higher temperatures was not measured. When polymerase activity was measured at temperatures <60 °C, the activity decreased by a mean of 6.7% every 5°C to 35°C. That being said, at temperatures 25 and 15°C the LR Bst is not as effective. However, the positive remains for it to be used in NA amplification protocols optimized for such temperatures.

Figure 1. Normalised Invitrogen Lyo-ready Bst DNA Polymerase activity at different temperatures



Invitrogen Lyo-ready Bst DNA Polymerase (LR Bst) activity testing at different incubation temperatures normalized to the activity of the LR Bst at 60°C. 3 replicates per test, results represent mean ± SD.


Overall, the LR Bst performed best at higher temperatures further demonstrating its suitability for reactions that require thermostable strand-displacement polymerases. LR Bst also maintained about half of its activity at 35-60°C which, considering its effective use in (RT-)LAMP at 3 U/rtxn, shows promise in low temperature reactions. Finally, some activity was still observed at <30°C that, with some optimisation, could be applied for diagnostics at room temperature.

Results

Invitrogen Lyo-ready Bst DNA Polymerase in LAMP: Speed

Bst polymerases are most commonly used in LAMP reactions, thus LAMP was chosen as the model NAAT to evaluate Bst polymerase functionality. Reaction speed was tested before any inhibitor testing was carried out. In all cases the Invitrogen Lyo-ready Bst DNA Polymerase exhibited the fastest reaction speed, allowing to go forward with inhibitor testing in LAMP (Figure 2).

Figure 2. Target detection time in LAMP and RT-LAMP




Reaction LAMP (A) and RT-LAMP (B) results for 30 min reactions at 65°C. 4 replicates per reaction. Dotted bars 475% variability. Error bars based on relative standard deviation (RSD). Specificity of results on heat control. A) DNA targets, M. pneumoniae using 4 primers only, no detection over 30 min = no time to detection shown. B) RNA targets with suggested 50.

LAMP in the presence of inhibitors

Loop-mediated isothermal amplification is one of the most popular current uses for Bst DNA polymerases in diagnostics. However, samples and their preparation can introduce a variety of inhibitors into a reaction. Thus, it is crucial that the LAMP reaction withstands a considerable amount of inhibiting material and still provides reliable results. Here, we selected a number of inhibitors that may be introduced through different means: ethanol (sample prep), KCl (sample or sample prep), humic acid, urea (sample) (Figure 3, 4, 5). The concentrations of inhibitors were chosen by probability to be introduced into the reaction.

Figure 3. LAMP results for samples containing ethanol.



Reaction LAMP: Target: 1000 copies Adv141, +4, 1% and 7% represent the total concentration of ethanol. Bars with no time to detection represent no detection within 30 min.

Results

Figure 4. LAMP results for samples containing higher KCl concentrations.





Figure 5. LAMP results for samples containing urea and humic acid.



Real-time LAMP: Target: 1000 copies Adv141, +4, 1% and 7% represent the total concentration of ethanol. Bars with no time to detection represent no detection within 30 min.

Overall, the Invitrogen Lyo-ready Bst DNA Polymerase, when compared to other commercially available Bst DNA polymerases, had faster target detection speed both with and without inhibitors in the reactions. Most importantly, Invitrogen Lyo-ready Bst DNA polymerase showed great resistance to salt concentration even at 175 mM additional KCl. High salt resistance can be correlated to higher processivity¹⁴ and help in reactions where high salt concentrations are present.

Conclusions

Invitrogen Lyo-ready Bst DNA Polymerase shows great promise in future diagnostic applications. Its wide active temperature range, salt tolerance, along with the exceptional functionality in LAMP reactions and lyophilisation compatible formulation allows for a variety of diagnostic uses and applications.

References

- Rahman MT, Sobur MA, Islam MS, Levy S, Hossain MJ, El Zowalaty ME, Rahman AT, Ashour HM. Zoonotic Diseases: Etiology, Impact, and Control. *Microorganisms*. 2020 Sep 12;8(9):1405. doi: 10.3390/micro8091405.
- Carlson CJ, Farrell MJ, Grange Z, Han BA, Molentzen N, Pheasant AL, Ramussen AL, Abery GF, Bett B, Brett-Major DM, Cohen LE, Dallas T, Eskew EA, Fague AC, Forbes KM, Gibb R, Hatalsi S, Hamner CC, Katz R, Kintrachuk J, Moyaert RL, Nutter FB, Ogata J, Oviak KJ, Rouke M, Ryan SJ, Ross N, Seifert SN, Sironen T, Stanley CJ, Taylor K, Verter M, Wehabe PW. The future of zoonotic risk prediction. *Philos Trans R Soc Lond B Biol Sci*. 2021 Nov 8;376(1837):20200358. doi: 10.1098/rstb.2020.0358. Epub 2021 Sep 20.
- Lindahl JF, Grace D. The consequences of human actions on risks for infectious diseases: a review. *Infect Ecol Epidemiol*. 2015 Nov 27;5:30048. doi: 10.3402/ieo.v5i30048.
- Kwok HF. The significance of advanced COVID-19 diagnostic testing in pandemic control. *Int J Biol Sci*. 2022 Jul 11;18(12):4610-4617. doi: 10.7150/ijb.72837.
- Augustine R, Hasan A, Das S, Ahmed R, Mori Y, Notomi T, Kevadiya BD, Thakor AS. Loop-Mediated Isothermal Amplification (LAMP): A Rapid, Sensitive, Specific, and Cost-Effective Point-of-Care Test for Coronavirus in the Context of COVID-19 Pandemic. *Biology (Basel)*. 2020 Jul 22;9(1):192. doi: 10.3390/biology909192.
- Administration, U. F. A. D. (2023). *Molecular Diagnostic Tests for SARS-CoV-2*. [online] Available at: <https://www.fda.gov/medical-devices/covid-19-emergency-use-authorizations-program/iv-diagnostic-assays-molecular-diagnostic-tests-sars-cov-2> [accessed 06/09/2023].
- Song X, Coulter, F.J., Yang, M. et al. A lyophilized colorimetric RT-LAMP test kit for rapid, low-cost, at-home molecular testing of SARS-CoV-2 and other pathogens. *Sci Rep* 12, 7043 (2022). <https://doi.org/10.1038/s41598-022-11144-5>
- Ziroz, P.G., Kokkinos, P.A., Alabd, A. et al. Development and Evaluation of a Loop-Mediated Isothermal Amplification Assay for the Detection of Adenovirus 40 and 41. *Food Environ Virol* 7, 279-286 (2015). <https://doi.org/10.1007/s12560-015-9182-8>
- Arfaatabar, M., Noori Goodarzi, N., Afshar, D., Memariani, H., Azimi, G., Masoorian, E. And POURMAND, M., 2019. Rapid Detection of Mycoplasma pneumoniae by Loop-Mediated Isothermal Amplification (LAMP) in Clinical Respiratory Specimens. *Tehran Journal of Public Health*.
- Fujino M, Yoshida N, Yamaguchi S, Hosaka N, Ota Y, Notomi T, Nakayama T. A simple method for the detection of measles virus genome by loop-mediated isothermal amplification (LAMP). *J Med Virol*. 2005 Jul;76(3):406-13. doi: 10.1093/jmv/2005/07.
- Alekseenko A, Barrett D, Parise-Sanchez Y, Howard RJ, Stradbruch E, Ampah-Korah H, Revanki N, Zuniga-Velaz S, Klenov A, Mallov A, Ye S, Liu X, Reiniss B, Elsäßer SJ, Nyman T, Sandt G, Yin X, Plechano V. Direct detection of SARS-CoV-2 using non-commercial RT-LAMP reagents on heat-inactivated samples. *Sci Rep*. 2021 Jan 18;11(1):1920. doi: 10.1038/s41598-020-80362-4.
- Lage JM, Leamon JH, Pejovic T, Hammons S, Lacey M, Dillon D, Segrares R, Vosbrinck B, Gonzalez A, Pinkel D, Albertson DG, Costa J, Lizardi PM. Whole genome analysis of genetic alterations in small DNA samples using hyperbranched strand displacement amplification and array-CGH. *Genome Res*. 2003 Feb;13(2):294-307. doi: 10.1101/jgv.377203.
- Beats, T.P., Smith, J.H., Ntupski, R.M. et al. A mechanism for amplified rolling circle amplification. *BMC Molecular Biol* 11, 94 (2010). <https://doi.org/10.1186/1471-2199-11-94>
- Notomi T, Okayama H, Masubuchi H, Yokoyama T, Watanabe H, Amino N, Hase T. Loop-mediated isothermal amplification of DNA. *Nucleic Acids Res*. 2000 Jun 16;28(12):E53. doi: 10.1093/nar/28.12.E53.
- Zhai B, Chow J, Cheng G. Two Approaches to Enhance the Processivity and Salt Tolerance of Staphylococcus aureus DNA Polymerase. *Protein J*. 2019 Apr;38(2):190-198. doi: 10.1007/s10930-019-09818-7.

Acknowledgements

We would like to thank Aistė Polikaitė, Valerija Prokopenė and Ingrida Vėnelė for valuable insights on isothermal amplification reactions during early research.

Trademarks/Licensing

© 2019 ThermoFisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.